

## Sanmedi Stoma Washbasin

Looking after a stoma requires quite a lot of attention. To make things easier for the stoma user, a comfortable washbasin has been designed. This washbasin is an excellent aid for stoma users and offers optimal hygiene.

### Comfortable

An individually adaptable stoma washbasin makes looking after the stoma easier and offers optimal hygiene. With the help of the stoma washbasin the user can do all the things he or she needs to do, like emptying the bag, changing the adhesive and rinsing the stoma in a comfortable and save way.

### For whom?

The stoma washbasin has been specially designed for colostomy, ileostomy and urostomy users. Especially for stoma users for whom irrigation is important the washbasin offers a modern and entirely new nursing method.

In case of individual installation the height of the dressing table is adapted to the height of the user which means that the upper side of the table is approximately 10 cm below the stoma. For general areas like in hospitals and institutions the standard height is between 80 and 85 cm.



### Ergonomic shape

The stoma washbasin has an ergonomic design. The front side is carefully adapted to the user's body, this means it is slightly lowered and curved to match the body. The front has a narrow edge which makes it very easy to hang the pouch inside the washbasin.

The stoma washbasin is made of very smooth, white ABS with a top coating. This material is warm to the touch. The one-handle mixing tap has a detachable hand-shower for extra cleaning purposes. The mirror is slanting so that it is easy to check the stoma. The powerful lateral flushing system ensures that all impurities are flushed away immediately.

### Where?

The stoma dressing table can be installed everywhere near a water supply and drainage point. That means it can be installed in private homes, in offices and institutions employing stoma users but also in public buildings like stoma clinics, hospitals and nursing homes.

For more information:

Tel. +31 251 316482

Fax +31 251 314043

E-mail: [info@sanmedi.nl](mailto:info@sanmedi.nl)